

Spa Range

VPN Bedpan - Hospital style bedpan with lid

Long Life Reusable Polypropylene Bedpan for use in Hospitals, Nursing, Rehabilitation and Elderly Care



The Spa Range is a range of re-usable living aids designed for use in Hospital Wards, Nursing and Care Homes, Rehabilitation and Personal Use at home.

Our Spa Range is an attractive, stylish and functional group of products which includes drinking cups, beakers, jugs, vomit bowls, medicine cups, urinals, bedpans and commodes.



Product Benefits

- Suitable for both male and female use.
- Ideal for patients confined to bed.
- Lightweight, strong and durable.
- Secure fitting lid to prevent leaks and spills and reduces cross infection and unpleasant odours.
- Can be cleaned in automatic washer disinfectors or hand washed.
- Smooth material surfaces ensure ease of use and ease of cleaning.
- Suitable for bariatric use - weight tested up to 250kg (39.5 stone).
- Guaranteed for a minimum of 1000 wash cycles.

Product Specification

- Bedpan with Lid - 2500ml capacity
- 335mm x 308mm x 108mm
- Pack Size 15
- Pack Weight 11.1 kg
- The following standards are applicable:
 - **ISO 9001:2008**
Quality Management System
 - **ISO 13485:2003**
Medical Devices - Quality Management System
 - **BS2588: 1995**
Re-usable Portable Sanitary Pans
- Manufactured from Polypropylene
- Maximum temperature 135 °C

Guide to Using Spa Products from Warwick SASCo

Usage

- Suitable for male and female use.
- For use in a seated position.
- Close lid after use to reduce odour and risk of spills.
- Empty after each use.
- Weight endurance tested up to 250 kg (39.5 stone). Advisable maximum usage time of up to 5 minutes - care should be given to patient comfort and safe positioning.

Cleaning

- **For Hospital Use** - Wash before first use. Clean in automatic washer disinfectors compliant to ISO 15883.
- Wash cycle at 80 °C. Disinfection cycle at 90 °C. For further processing information, refer to machine manufacturers instructions and local hospital policy.
- **For Personal Use** - Wash before first use in hot water with detergent.
- Empty contents after use down a toilet. Rinse pan empty down a toilet.
- Wash in hot water with detergent after each use and allow to dry naturally.

Storage

- This product has a shelf life of 36 months from date of manufacture providing the recommended storage conditions are adhered to.
- Product should be stored in its primary packaging, off the floor on a pallet or shelf, and in a clean dry ventilated room at a temperature of 4 °C - 25 °C.