

SteriLite Range

GPS120 Gallipot

Single Use Disposable Polypropylene Gallipot for use in Disposable Sets, Decontamination Services, Hospitals and Clinics.



The SteriLite Range comprises of single use, disposable products and was introduced by Warwick SASCo to speed up packing and sterilizing procedures, reducing the need for retrieval and decontamination.

It is ideal for use in Sterile Services, Decontamination Units, Operating Rooms and Clinics due to its low cost, product design, and performance.

The product range now offers single use, disposable trays, kidney dishes, gallipots, sponge bowls, quivers, jugs and medicine measures.



Product Benefits

- Lightweight, strong and durable.
- Ideal for use in single use surgical and custom packs and hospital processing.
- Clear graduations for easy measuring.
- Low cost.
- Clean and ready for sterilization.
- Can be sterilized using Autoclave, Ethylene Oxide or Hydrogen Plasma.
- Easy release.
- Smooth edges and surfaces eliminate paper rips during packing and sterilization
- Autoclavable up to 135°C
- Non-sterile but clean and easy for packing and sterilizing

Product Specification

- 120ml Gallipot
- 85mm x 30mm
- Pack Size 1000
- The following standards are applicable:
 - **ISO 9001:2008**
Quality Management System
 - **ISO 13485:2003**
Medical Devices - Quality Management System
- Manufactured from high quality Polypropylene
- SteriLite products are single use and should not be re-processed after use.
- SteriLite products are provided non-sterile but clean and ready for packing and sterilization.
- Keep SteriLite products stored in their polythene bags to avoid contamination.
- This product has a self life of 36 months from date of manufacture providing the recommended storage conditions are adhered to.

Guide to Using SteriLite Products from Warwick SASCo

A standard packing material may be used. Ensure that the pack is large enough to contain the device without stressing seals. We suggest the following:

- A light weight pouch material (60gms)
- Pouch to be big enough to allow shrinkage round product choose a pouch where the size of the device does not exceed 75% of the space inside the pouch.
- Choose a pouch with a gusset (50mm at least)
- The rim of the GPS120 to be facing the white paper
- Allow plenty of space between product and pleat in pouch
- Do not force the GPS120 into the pouch ensure there is plenty of space between products
- Single wrapping is preferable
- Ensure that the pack is large enough to contain the device without stressing seals
- Product should be stored in its primary and secondary packaging, off the floor on a pallet or shelf, and in a clean dry ventilated room at a temperature of 4 °C - 25 °C.