

SURGICAL INSTRUMENT PURCHASE AND CARE GUIDE

Guidance for purchasing and caring for surgical instruments
from ABHI's Surgical Instruments Special Interest Section Group



CONTENTS

The Purpose Of This Guide	2
Understanding Quality	
Buying The Right Instrument Is A Collective Responsibility	3
Know Your Standards	4
Know Your Materials	5
Manufacturing Processes	6
Common Features And Terminology For Quality Instruments	7
Identifying Common Instrument Features	8-9
Maintaining High Quality	
A Guide To Reprocessing Re-Usable Surgical Instruments	10-11
Care And Maintenance Tips	12-13
Ethical Supply	14

THE PURPOSE OF THIS GUIDE

This booklet is designed to help healthcare providers achieve the best whole life value for money in their purchasing decisions.

Surgical instruments are a critical component of surgical procedures. It is important that purchasers are well informed, to ensure patient safety as well as best value.

This Guide is an educational and training tool. It helps improve awareness and understanding of how surgical instruments are made, the standards which apply to them and the quality of the instruments.

By enabling effective procurement, we hope to help healthcare providers achieve the best return on their investment, while putting patients at the heart of decision-making.

The Association of British Healthcare Industries (ABHI) is the UK's industry association for the medical technology sector. The companies we represent produce around 85% of the industry's total UK output.

We promote the rapid adoption of medical technologies in the UK and key global markets to maximise patient outcomes, and support ethical procurement.


“ Quality is always top and non-negotiable ”

Lord Carter Health and Care Show, July 2016

UNDERSTANDING QUALITY

BUYING THE RIGHT INSTRUMENT IS A COLLECTIVE RESPONSIBILITY

Purchasing Surgical Instruments needs to be a co-ordinated process with input from the appropriate health professionals before and after purchase:



The surgeon, theatre staff, sterilisation and decontamination teams are all essential to surgical instrument purchasing decisions. Their feedback is critical in making the right decision.

Health Technical Memorandum 01-01: Management and decontamination of surgical instruments (medical devices) used in acute care (or as appropriate)

UNDERSTANDING QUALITY

KNOW YOUR STANDARDS

Surgical instruments are governed by a number of standards including, but not limited to:

Medical Devices Directive 93/42/EEC

- MDD – this Directive includes the essential requirements such as CE marks to be followed by manufacturers.

ISO 7153-1:2001 BS 51994-1:1991

- The Standard for the composition of the different materials and steel grades used.

BS 5194-4:1985

- For the specifications of scissors, shears, and other cutting instruments.

ISO 13485

- Requirements for a quality management system, where an organisation needs to demonstrate its ability to provide medical devices.

BS 5194-4:1989

- For the specifications of instruments with pivot points.

BS 5194-3:1995

- For the specifications of dissecting forceps.

CE Marking

- On every device, look for a CE mark, the name of the manufacturer and a traceability code. Be aware that a CE mark is a sign of compliance with MDD however and should not be taken as an automatic sign of quality.



UNDERSTANDING QUALITY

KNOW YOUR MATERIALS

Most surgical instruments start life as forgings or “blanks”. They are governed by two International Standards for material specification: DIN 17442 and DIN EN 10088-3 8/95.

Surgical instruments are mainly made from two types of stainless steel: martensitic and austenitic. Some are made from titanium. The boxes on the right illustrate the types of instrument materials.

- ▶ **Martensitic is magnetic and contains up to 1% carbon which allows the instrument to be heat-treated**
- ▶ **Austenitic is the most common type of stainless steel and is highly versatile**

ISO 7153-1 has a full list of the suitable grades of stainless steel available.

Martensitic Grade B-420 S29

- ▶ **Used for non-cutting instruments, e.g. artery forceps**

Hardness 40-48 HRC
Carbon content 0.16-0.25%
Chromium content 12-14%

Martensitic Grade C or D-420 S45

- ▶ **Used for cutting instruments, e.g. scissors & gouges**

Hardness 50-58 HRC
Carbon content 0.35-0.45%
Chromium content 12-14%

Austenitic Grade 304 S15

- ▶ **Used for instruments which do not require hardening, e.g. dental tweezers and holloware**


Hardness 40-48 HRC
Carbon content 0.07-0.15%
Chromium content 16-19%
Nickel content 8-11%

Titanium

- ▶ **Used for Ophthalmic & Microsurgery instruments**

Ti-6Al-4V ELI or grade 23 titanium


MANUFACTURING PROCESSES



COMMON FEATURES AND TERMINOLOGY FOR QUALITY INSTRUMENTS

There are a huge variety of features which appear on reusable surgical instruments. Here are a few common features and what to look for in a quality item:

- ▶ Teeth and prongs should be sharp and mesh exactly when jaws close
- ▶ Serrations on both jaws should be identically shaped and mesh exactly
- ▶ When pressure is released, the teeth, prongs, and serrations should part freely without catching
- ▶ No slippage in Needleholder jaws
- ▶ Should be symmetrical
- ▶ Rectangular section should give maximum strength to the joint
- ▶ Should avoid unnecessary gaps
- ▶ Use of countersink prevents rivet from moving
- ▶ Joint should move smoothly, not too tight, not too loose
- ▶ It should be possible to open and close the joint easily with 2 fingers



IDENTIFYING COMMON INSTRUMENT FEATURES & WHERE THEY MAY BE SEEN

Box Joint

Found on:

- Spencer Wells artery forceps
- Halstead Mosquito artery forceps
- Crile artery forceps



Screw Joint

Found on:

- Mayo scissors
- Metzenbaum scissors
- McIndoe scissors
- Stevens scissors



Ratchet

Found on:

- Norfolk and Norwich retractors
- Travers retractors
- West and Weitlander retractors



Teeth

Found on:

- Littlewoods tissue forceps
- Allis tissue forceps
- Lanes tissue forceps



Rack

Found on:

- Mayo-Hegar needle holder
- Spencer Wells artery forceps



Scissor Blades

Found on:

- Mayo scissors
- Metzenbaum scissors
- Dressing scissors



Lumen

Found on:

- ▶ American, Adson, Magill & Zoellner suction tubes



Finish

- ▶ Bright polished or satin finished steel



Tungsten Carbide Tips

Found on:

- ▶ Needle holders for durability



Serrations

Found on:


- ▶ Artery forceps



Atraumatic Teeth

Found on:

- ▶ Debakey clamps and forceps
- ▶ Derra & Cooley vascular clamps




Bows

Found on:


- ▶ Scissors, needle holders & artery forceps



MAINTAINING HIGH QUALITY A GUIDE TO REPROCESSING RE-USABLE SURGICAL INSTRUMENTS



To prepare for decontamination, reprocess all instruments as soon as practicable following use. Disassemble only where intended, without the use of tools, unless specifically made available with the instrument and provided by the manufacturer.



Automated Cleaning

Use CE marked or validated washer-disinfector machines and low-foaming, non-ionic cleaning agents and detergents. Follow the manufacturers' instructions for use, warnings, concentrations, and recommended cycles.

- Load instruments carefully, with box joints and hinges open, and so that any fenestrations can drain.
- Place heavy instruments with care in the bottom of containers. Do not overload wash baskets.
- Place instruments with concave surfaces facing down to prevent pooling of water.
- Use appropriate attachments to flush in side reamers, and devices with lumens or cannula.
- Ensure that soft, high purity water which is controlled for bacterial endotoxins is used in the final rinse stage.

Note Automated Cleaning may not be suitable for all lumens and cannula, in which case clean manually with a water jet gun, if available, and an appropriate brush and/or stilette that reaches the depth of the feature.

After manually cleaning, pass all devices through an automatic cleaning cycle to achieve disinfection.

FROM POINT OF USE

PREPARATION

CLEANING

Wherever possible, do not allow blood, debris or bodily fluids to dry on instruments. To prolong their life, reprocess immediately after use. If that's not possible, use an enzymatic foam spray to help prevent soil from drying.



IMPORTANT NOTE: This is not comprehensive. For a full, validated reprocessing guide, speak to your instrument supplier and follow current MHRA guidelines for reprocessing instruments.



INSPECTION

After cleaning, visually inspect all



- ▶ Surfaces
- ▶ Ratchets
- ▶ Cannulations
- ▶ Holes
- ▶ Joints

Lumens for complete removal of soil and fluids. If any soil or fluid is still visible, return the instrument for repeat decontamination.

Ensure that instruments are dry before sterilisation. Always follow the instructions of the machine manufacturer. Use a CE marked or validated vacuum autoclave operating at 134-137°C 2.25 bar for 3 minutes minimum holding time.

When sterilising multiple instruments in one cycle, always make sure that the stated maximum load is not exceeded.

PACKAGING & STERILISATION



STORAGE

All instruments to be packed following local protocol in accordance with BS standards.

MAINTAINING HIGH QUALITY CARE & MAINTENANCE TIPS

For full guidance see www.a-k-i.org 'Red Brochure'

BROKEN/CRACKED BOX JOINTS

Tension stress

Cause Heating and cooling in sterilisation process

Solution Close instrument to first notch only during sterilisation



Forced stress

Cause Overloading instruments

Solution Ensure correct device and attachment is being used

NB: Also sutures and needle holders



General stress

Cause Build up of blood and debris in box joint

Solution Ensure instruments are cleaned in open position during washing and disinfection



DISCOLOURATION

Water spots

Light coloured, often with sharply defined edges

Cause Final rinse or sterilisation water supply contains high concentration of minerals

Solution Use demineralised water in final rinse, and pure steam in sterilisation



Oxidisation spots

Light coloured, without sharply defined edges

Cause Final rinse or sterilisation water supply contains high concentration of heavy metal ions and/or silicates

Solution Use demineralised water in final rinse, and pure steam in sterilisation



Yellow brown to dark brown spots

Cause Debris has dried on the device before cleaning or hasn't been removed due to poorly performing detergents

Solution Remove by thoroughly scrubbing with a good detergent, otherwise corrosive pitting will occur




Other causes of discolouration

- ✦ Insufficient rinsing off detergents and disinfectants
- ✦ Chlorides
- ✦ Water droplets slowly condensing on instruments during sterilisation
- ✦ Inferior detergent

CORROSION

Pitting corrosion

Cause 1 Excessive chloride concentrations
Solution Use demineralised water



Cause 2 Prolonged exposure to saline solutions (blood, debris or contaminated disinfectant or detergent) where bacterial activity creates acidic residue
Solution Clean instruments as soon as possible after use

Abrasion corrosion


Cause Build up of debris stops devices from opening and operating smoothly, causing destruction of passivation layer at joints and crevices



Solution Ensure instruments are cleaned in open position & lubricate regularly

Contaminated steam corrosion

Cause Rusty steam in sterilisation process
Solution Regular validation and maintenance of decontamination equipment



Surface corrosion

Cause Damage to passivation layer
Solution Avoid use of strong acid, alkaline or caustic solutions
NB: Aluminium is particularly susceptible



Spreading corrosion

Cause Instruments sterilised with already rusty devices - rust is transferred through the detergent solutions
Solution Separate rusty devices from "healthy" ones



The Red Brochure

The Instrument Reprocessing Working Group was set up in 1976. They have produced a Surgical Instrument guidance document for the past 40 years. This provides exhaustive guidance on all aspects of surgical instrument care and best practice. www.a-k-i.org



ETHICAL SUPPLY

Ethical Manufacturing & the NHS Supply Chain's Labour Standards Assurance System (LSAS)

ABHI has its own code of business practice and we support the ethical sourcing of products. The Surgical Instruments SIS Group worked with NHS Supply Chain as part of the 2012 (and pending 2017) Surgical Instruments Framework Agreement to launch its Labour Standards Assurance System.

LSAS is a matrix of ethical requirements designed by NHS Supply Chain and the Department of Health, through which suppliers are audited and assessed by a third party notified body. The responsibility is with the supplier to ensure there is continual progress and regular risk assessment and review, to mitigate potential ethical and labour risks in the supply chain.

This has been embedded since 2012 and many of our members have improved to obtain level 2 and 3 on the framework.

ABHI is committed to promoting good ethical practice amongst members, we see this as integral and essential for improving labour standards in both single use and reusable surgery instrument manufacturing.

ABHI Code of Business Practice

At ABHI, we place ethical compliance at the heart of the medical technology industry. Healthcare professionals and patients must feel they can be confident in our ethical standards at all times, so they can work with us to improve the innovations we develop.

We have been working hard for several years to help member and other companies reach the highest standards – both as organisations and as individuals at all levels.

It is a condition of ABHI membership that a company adheres to the ethical standards in the ABHI Code of Business Practice. The Code stipulates minimum standards for members' business practices in the UK, Europe and elsewhere.

✦ **More information can be found at**
www.abhicodeofpractice.org.uk





**This document has been produced
by members of the ABHI Surgical
Instruments Special Interest Section.**

ABHI Member details

Association of British Healthcare Industries
107 Gray's Inn Rd, London WC1X 8TZ
Tel: +44 (0)20 7960 4360 Fax: +44 (0)20 7960 4361
enquiries@abhi.org.uk www.abhi.org.uk

 @UK_ABHI

Version 1

