

# PRODUCT TECHNICAL DATASHEET

## Product Profile & Instructions for Use

Warwick House, Heathcote Way, Warwick  
CV34 6TE, United Kingdom  
T: 0044 (0)1926 422427 [www.sasco.co.uk](http://www.sasco.co.uk)  
[enquiries@sasco.co.uk](mailto:enquiries@sasco.co.uk)

### TR30 Attachment Ring

**Not a medical device**

**Dimension : 30 mm**

**Material : Stainless Steel**

**Pack Size : 1 x 100**



#### Intended Use

- Attachment ring for use in sterile services and decontamination, for identification and traceability.
- To attach an identification tag to a tray.
- Reusable product.

#### Intended User

- Product to be used by clinical professionals and decontamination technicians.

#### Product Features

- Autoclavable.
- Can securely attach tags to instrument trays and containers to track and trace through the decontamination process.
- Guaranteed 1000 cycles for steam sterilization.
- Aids clear and easy identification of sets during processing.
- Lightweight, strong, and durable.
- Smooth edges and surfaces eliminate paper rips during sterilization.

#### Packaging Information

Catalogue Number	Case Bar Code (GTIN 14)	Primary Packaging	Secondary Packaging	Pack Weight
TR30	15060098131538	Polythene bag 1 x 100	n/a	0.25 Kg

#### Product Symbols

Manufacturer	Catalogue number	Caution	Batch code	Use-by date	Non-sterile	Absence of natural rubber latex	Do not use if package is damaged	Keep dry	Keep away from sunlight

#### Additional Information

- Harmonised System – Product Tariff / Commodity Export Code 7326 2000 90.
- Country of Origin United Kingdom.

#### Manufacturer

- Warwick SASCo Ltd. Warwick House, Heathcote Way, CV34 6TE, United Kingdom.
- Telephone: 0044 (0)1926 422427 E-mail: [enquiries@sasco.co.uk](mailto:enquiries@sasco.co.uk)



## Product Profile & Instructions for Use

Warwick House, Heathcote Way, Warwick  
CV34 6TE, United Kingdom  
T: 0044 (0)1926 422427 www.sasco.co.uk  
enquiries@sasco.co.uk

### Distributor/Importer



### Standards

- **Company compliant to:** **BS EN ISO 13485:2016 Medical Devices. Quality management systems. Requirements for regulatory purposes.**  
**BS EN ISO 14971:2019 Medical Devices. Application of risk management to medical devices.**  
**BS EN ISO 15223-1:2021 Medical Devices. Symbols to be used with information to be supplied by the manufacturer. Part 1: General requirements.**  
**BS ISO 17664-2:2021 Processing of health care products – Information to be provided by the medical device manufacturer for the processing of medical devices. Part 2: Non-critical medical devices.**  
**BS EN ISO 20417:2021 Medical devices. Information to be supplied by the manufacturer.**
- **Referenced standards:** **BS EN ISO 15883 Washer disinfectors. Requirements.**  
**BS EN ISO 11135 Sterilization of health care products. Ethylene oxide.**  
**BS EN ISO 17665-1 Sterilization of health care products. Moist heat.**  
**BS ISO 22441 Sterilization of health care products. Low temperature vaporized hydrogen peroxide.**

## Instructions for Use

### Storage and Handling

- **Do not drop.**
- **Store in original packaging.**
- **Store in clean and dry conditions.**
- **Store off the floor on a shelf or pallet.**
- **Shelf life 36 months.**

### Preparation

- **Product supplied non-sterile.**
- **Inspect products before processing.**
- **Withdraw damaged products and report to manager.**
- **Rinse with warm water to remove any packaging residues.**

### Product can be used in conjunction with

- **Tray Tag; TAG (all colours).**
- **Instrument Tray – Mesh base; MIT1515, MIT2709, MIT2715, MIT2817, MIT3025, MIT4030**
- **DIN Style Tray – Mesh base and sides; DINT2525, DINT4825**

### Cleaning and Disinfection

- **Position products individually, avoiding contact with other items.**
- **Position products to allow water penetration and drainage.**
- **Machine clean in washer disinfectors compliant with ISO 15883.**



## Product Profile & Instructions for Use

Warwick House, Heathcote Way, Warwick  
CV34 6TE, United Kingdom  
T: 0044 (0)1926 422427 www.sasco.co.uk  
enquiries@sasco.co.uk

- Use approved detergents, enzymatic or biocide cleaners, suitable for cleaning medical devices and surgical instruments which are compatible with stainless steel material (see detergent manufacturers' guide).
- Use detergents and disinfectants with a range between 2pH and 13pH.
- Ultrasonic cleaning is permissible.
- Manual cleaning is permissible – separate product components and clean individually.
- Rinse with clean water.
- Dry with unwoven towels.
- Do not use abrasive materials.

### Inspection and Maintenance

- Visually inspect product to identify damage.
- Discard damaged product.

### Wrapping

- Sterilization wrapping sheets and pouches must be compatible with stainless steel material and suitable for the chosen method of sterilization.
- Be aware that some wrapping sheets and pouches shrink during sterilization.
- Pouches must be large enough to contain the product without stressing the seals.
- We advise using a pouch or bag that is 25% larger than the product.

### Sterilization

- Moist Heat sterilization (Autoclave) at 134.6 °C (+2 °C) for a minimum of 2 minutes.
- EO Ethylene Oxide "Gas" sterilization at maximum temperature of 65 °C.
- H2O2 Hydrogen Peroxide "Plasma" sterilization at maximum temperature of 65 °C.
- Sterilization procedures and machines to comply with ISO 17665-1, ISO 22441, and ISO 11135.

### Disposal

- End of life is determined by wear and tear, visual inspection, or damage to product.
- If product is contaminated treat as clinical waste, and dispose of in accordance with clinical waste procedures.

### Warnings

- Do not use abrasive powders or metal brushes as they may cause surface damage.
- Do not exceed 137 °C for steam sterilization.
- Do not exceed 8 minutes hold time for steam sterilization.

### Important Information for Processing

- The information provided by Warwick SASCo is the company's official instructions for this product.
- It is the responsibility of the processor to ensure the processing is performed in accordance with this information and that personnel, machines and equipment comply with relevant standards for processing.
- Any deviation from this information requires evaluation for the effectiveness and potential adverse consequences.
- Additional information available on request contact: enquiries@sasco.co.uk

