

LCT2215S Packing Tray - 4 compartments blue

Single Use Disposable Polypropylene Compartment Tray for use in Disposable Sets, Decontamination Services, Hospitals and Clinics.



The Sterilite Range comprises of single use, disposable products and was introduced by Warwick SASCo to speed up packing and sterilizing procedures, reducing the need for retrieval and decontamination.

It is ideal for use in Sterile Services, Decontamination Units, Operating Rooms and Clinics due to its low cost, product design, and performance.

The product range now offers single use, disposable trays, kidney dishes, gallipots, sponge bowls, quivers, jugs and medicine measures.



Warwick House, Heathcote Way, Warwick, CV34 6TE

T: +00 44 (0)1926 422427

F: +00 44 (0) 1926 315405

E: info@sasco.co.uk

W: www.sasco.co.uk

Product Benefits

- Lightweight, strong and durable compartment tray with integral gallipots to maximise space for instruments and accessories.
- Ideal for use in single use surgical and custom packs and hospital processing.
- Low cost.
- Can be sterilized using Autoclave, ETO or Hydrogen Plasma.
- Grip feature moulded below the rim to improve holding.
- Smooth edges and surfaces eliminate paper rips during packing and sterilization.

Product Specification

- 220mm x 150mm x 24mm
- Pack Size 160
- The following standards are applicable:
 - **ISO 9001:2008**
Quality Management System
 - **ISO 13485:2003**
Medical Devices - Quality Management System
- Manufactured from high quality Polypropylene
- Non-sterile but clean and ready for packing
- Autoclavable up to 135 ° C

Guide to Using Sterilite Products from Warwick SASCo

- Sterilite products are provided non-sterile but clean and ready for packing and sterilization.
- Keep Sterilite products stored in their polythene bags to avoid contamination.
- Sterilite products are single use and should not be re-processed after use.
- This product has a shelf life of 36 months from date of manufacture providing the recommended storage conditions are adhered to.
- Product should be stored in its primary and secondary packaging, off the floor on a pallet or shelf, and in a clean dry ventilated room at a temperature of 4 °C - 25 °C.