# Spa Range



## VB300 Vomit Bowl - With handle

Long Life Re-usable Polypropylene Vomit Bowl for use in Nursing, Rehabilitation, Elderly Care, and Community Care



The Spa Range is a range of re-usable living aids designed for use in Hospital Wards, Nursing and Care Homes, Rehabilitation and Personal Use at home

Our Spa Range is an attractive, stylish and functional group of products which includes drinking cups, beakers, jugs, vomit bowls, medicine cups, urinals, bedpans and commodes.

#### **Product Benefits**

- Ideal for patients confined to bed due to illness or surgery.
- Lightweight, strong and durable.
- Can be cleaned in automatic washer disinfectors or hand washed.
- Smooth material surfaces ensure ease of use and ease of cleaning.

## **Product Specification**

- Vomit Bowl with Handle 1.5 litre capacity
- 300mm x 144mm
- Pack size 25
- The following standards are applicable:
  - ISO 9001:2008

Quality Management System

- ISO 13485:2003

Medical Devices - Quality Management System

BS5452

Hospital Hollow-ware made of Plastics Material

Manufactured from high quality Polypropylene

## **Guide to Using Spa Products from Warwick SASCo**

#### Usage

- For use with nursing ill patients ideal for managing the after effects of anaesthesia.
- Empty after each use.

#### Cleaning

- Product should be washed before first use. Wash in hot water with detergents and disinfectants and allow to dry naturally.
- Wash immediately after in hot water with detergents and disinfectants.
- Spa Products are designed for repeated use guaranteed for 1000 washes.
- Scouring and abrasive cleaners will reduce the performance of products.

### Storage

- This product has a shelf life of 36 months from date of manufacture providing the recommended storage conditions are adhered to.
- Product should be stored in its primary packaging, off the floor on a pallet or shelf, and in a clean dry ventilated room at a temperature of 4 °C - 25 °C.

