

Spa Range

JGL1000 Jug with Lid - 1000ml

Long Life Re-usable Polypropylene Beaker for use in Nursing, Rehabilitation, Elderly Care and Community Care



The Spa Range is a range of re-usable living aids designed for use in Hospital Wards, Nursing and Care Homes, Rehabilitation and Personal Use at home.

Our Spa Range is an attractive, stylish and functional group of products which includes drinking cups, beakers, jugs, vomit bowls, medicine cups, urinals, bedpans and commodes.



Product Benefits

- Clear, easy to read graduated measurements.
- Lightweight, strong and durable.
- Easy to handle
- Ideal for use in Hospitals, Clinics, Nursing and Care Homes.
- Can be sterilized by Autoclave, Ethylene Oxide or Hydrogen Plasma.
- Can also be washed in automatic washer disinfectant or by hand using detergents and disinfectants.
- Also available with lid - JG1000

Product Specification

- 1000ml Measuring Jug with lid- Graduated Scale
- Available in Natural
- 142mm high
- Pack Size 40
- The following standards are applicable:
 - **ISO 9001:2008**
Quality Management System
 - **ISO 13485:2003**
Medical Devices - Quality Management System
- Manufactured from high quality Polypropylene

Guide to Using Spa Products from Warwick SASCo

Usage

- For use with hot and cold liquids.
- Do not overfill with liquids.
- Place lid firmly onto beaker, ensuring the spout centre is central to the beaker.

Cleaning

- Product should be washed before first use.
- Wash in hot water with detergents and allow to dry naturally.
- Spa Products are designed for repeated use - guaranteed for 1000 washes.
- Scouring and abrasive cleaners will reduce the performance of products.

Storage

- This product has a shelf life of 36 months from date of manufacture providing the recommended storage conditions are adhered to.
- Product should be stored in its primary packaging, off the floor on a pallet or shelf, and in a clean dry ventilated room at a temperature of 4 °C - 25 °C.