

Spa Range

TBC Beaker - Clear glass style

Long Life Re-usable Drinking Beaker for use in Hospitals, Nursing, Rehabilitation, Elderly Care, and Community Care



The Spa Range is a range of re-usable living aids designed for use in Hospital Wards, Nursing and Care Homes, Rehabilitation and Personal Use at home.

Our Spa Range is an attractive, stylish and functional group of products which includes drinking cups, beakers, jugs, vomit bowls, medicine cups, urinals, bedpans and commodes.



Product Benefits

- Lightweight, strong and durable.
- Ideal as a beaker for bathroom or bedside use.
- Clear plastic allows easy viewing of contents.
- Can be used for hot or cold liquids.
- Also available as part of our Locker Set (LST) - beaker, jug, lid and tray.

Product Specification

- Transparent Glass Style Beaker - 280ml
- Pack Size 120
- The following standards are applicable:
 - **ISO 9001:2008**
Quality Management System
 - **ISO 13485:2003**
Medical Devices - Quality Management System
- Manufactured from SAN.

Guide to Using Spa Products from Warwick SASCo

Usage

- For use with hot and cold liquids.
- Do not overfill with liquids.

Cleaning

- Product should be washed before first use. Wash in hot water with detergents and disinfectants and allow to dry naturally. Max temperature 110 °C.
- After use wash in hot water using detergents and disinfectants.
- Spa Products are designed for repeated use - guaranteed for 1000 washes.
- Scouring and abrasive cleaners will reduce the performance of products.

Storage

- This product has a shelf life of 36 months from date of manufacture providing the recommended storage conditions are adhered to.
- Product should be stored in its primary packaging, off the floor on a pallet or shelf, and in a clean dry ventilated room at a temperature of 4 °C - 25 °C.