

## MM60 Medicine Measure - Graduated

Long Life Re-usable Medicine Measure for use in Decontamination Services, Hospitals, Clinics and Care Homes



The Spa Range is a range of re-usable living aids designed for use in Hospital Wards, Nursing and Care Homes, Rehabilitation and Personal Use at home.

Our Spa Range is an attractive, stylish and functional group of products which includes drinking cups, beakers, jugs, vomit bowls, medicine cups, urinals, bedpans and commodes.



### Product Benefits

- Lightweight, strong and durable.
- Clear graduations for accurate measurement.
- For dispensing medication in hospitals, clinics and nursing homes.
- Can be used in custom packs.
- Can be sterilized by Autoclave, Ethylene Oxide or Hydrogen Plasma.

### Product Specification

- Medicine Measure graduated to 60ml
- Pack Size 1500
- Pack Weight 10kg
- Registered as Class 1 Medical Device (CE Marked)
- Complies with following Standards:
  - **ISO 9001: 2008**  
Quality Management System
  - **ISO 13485:2003**  
Medical Devices - Quality Management System
- Manufactured from polypropylene
- Autoclavable up to 135 ° C

### Guide to Using Warwick Products from Warwick SASCo

#### Usage

- Care should be taken when measuring liquid medicines as some medicines may stick to inside wall.
- Measuring tolerance +/- 2ml at 50ml or 60ml graduated.
- Critical dosages should be measured with a syringe.
- Place cup on level surface for accurate reading.

#### Cleaning

- Product should be washed before first use. Wash in hot water with detergent.
- Maximum processing temperature 135 ° C

#### Storage

- This product has a shelf life of 36 months from date of manufacture providing the recommended storage conditions are adhered to.
- Product should be stored in its primary packaging, off the floor on a pallet or shelf, and in a clean dry ventilated room at a temperature of 4 ° C - 25 ° C.